



NOTCH2 Antibody, HRP conjugated

Product Code	CSB-PA015950LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04721
Immunogen	Recombinant Human Neurogenic locus notch homolog protein 2 protein (59-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurogenic locus notch homolog protein 2 (Notch 2) (hN2) [Cleaved into: Notch 2 extracellular truncation (N2ECD); Notch 2 intracellular domain (N2ICD)], NOTCH2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NOTCH2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.