



# NTRK1 Antibody, FITC conjugated

|                     |  |
|---------------------|--|
| Product Code        | CSB-PA016133LC01HU   |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| Uniprot No.         | P04629   |
| Immunogen           | Recombinant Human High affinity nerve growth factor receptor protein (33-208AA)  |
| Raised In           | Rabbit   |
| Species Reactivity  | Human  |
| Form                | Liquid   |
| Conjugate           | FITC   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4   |
| Purification Method | >95%, Protein G purified   |
| Isotype             | IgG  |
| Clonality           | Polyclonal   |
| Alias               | High affinity nerve growth factor receptor (EC 2.7.10.1) (Neurotrophic tyrosine kinase receptor type 1) (TRK1-transforming tyrosine kinase protein) (Tropomyosin-related kinase A) (Tyrosine kinase receptor) (Tyrosine kinase receptor A) (Trk-A) (gp140trk) (p140-TrkA), NTRK1, MTC TRK TRKA |
| Immunogen Species   | Homo sapiens (Human)   |
| Research Area       | Others   |
| Target Names        | NTRK1  |