



NUMA1 Antibody, HRP conjugated

Product Code	CSB-PA016185YB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14980
Immunogen	Recombinant Human Nuclear mitotic apparatus protein 1 protein (627-870AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear mitotic apparatus protein 1 (Nuclear matrix protein-22) (NMP-22) (Nuclear mitotic apparatus protein) (NuMA protein) (SP-H antigen), NUMA1, NMP22 NUMA
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NUMA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.