



PEX1 Antibody, FITC conjugated

Product Code	CSB-PA017793LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43933
Immunogen	Recombinant Human Peroxisome biogenesis factor 1 protein (927-1283AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisome biogenesis factor 1 (Peroxin-1) (Peroxisome biogenesis disorder protein 1), PEX1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PEX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.