



# PEX1 Antibody, Biotin conjugated

Product Code	CSB-PA017793LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43933
Immunogen	Recombinant Human Peroxisome biogenesis factor 1 protein (927-1283AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisome biogenesis factor 1 (Peroxin-1) (Peroxisome biogenesis disorder protein 1), PEX1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PEX1