



PSD Antibody, HRP conjugated

Product Code	CSB-PA018842LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A5PKW4
Immunogen	Recombinant Human PH and SEC7 domain-containing protein 1 protein (247-489AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PH and SEC7 domain-containing protein 1 (Exchange factor for ADP-ribosylation factor guanine nucleotide factor 6) (Exchange factor for ARF6) (Exchange factor for ARF6 A) (Pleckstrin homology and SEC7 domain-containing protein 1), PSD, EFA6 EFA6A KIAA2011 PSD1 TYL
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PSD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.