



SERPINA5 Antibody, FITC conjugated

Product Code	CSB-PA021061LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05154
Immunogen	Recombinant Human Plasma serine protease inhibitor protein (157-367AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plasma serine protease inhibitor (Acrosomal serine protease inhibitor) (Plasminogen activator inhibitor 3) (PAI-3) (PAI3) (Protein C inhibitor) (PCI) (Serpin A5), SERPINA5, PCI PLANH3 PROCI
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SERPINA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.