



SV2A Antibody, Biotin conjugated

Product Code	CSB-PA022978LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L0J3
Immunogen	Recombinant Human Synaptic vesicle glycoprotein 2A protein (36-149AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Synaptic vesicle glycoprotein 2A, SV2A, KIAA0736
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SV2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.