





TBC1D4 Antibody

Product Code	CSB-PA023218LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60343
Immunogen	Recombinant Human TBC1 domain family member 4 protein (1-94AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TBC1 domain family member 4 (Akt substrate of 160 kDa) (AS160), TBC1D4, AS160 KIAA0603
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TBC1D4


Image



Immunohistochemistry analysis of human adrenal gland tissue using CSB-PA023218LA01HU at dilution of 1:100



Immunocytochemistry analysis of human kidney tissue using CSB-PA023218LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA023218LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.