



TBC1D4 Antibody, HRP conjugated

Product Code	CSB-PA023218LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60343
Immunogen	Recombinant Human TBC1 domain family member 4 protein (1-94AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TBC1 domain family member 4 (Akt substrate of 160 kDa) (AS160), TBC1D4, AS160 KIAA0603
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TBC1D4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.