



TBC1D4 Antibody, FITC conjugated

Product Code	CSB-PA023218LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60343
Immunogen	Recombinant Human TBC1 domain family member 4 protein (1-94AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TBC1 domain family member 4 (Akt substrate of 160 kDa) (AS160), TBC1D4, AS160 KIAA0603
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TBC1D4