




# TBX1 Antibody

|                            |                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA023245YA01HU                                                                      |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                           |
| <b>Uniprot No.</b>         | O43435                                                                                  |
| <b>Immunogen</b>           | Recombinant Human T-box transcription factor TBX1 protein (22-107AA)                    |
| <b>Raised In</b>           | Rabbit                                                                                  |
| <b>Species Reactivity</b>  | Human                                                                                   |
| <b>Tested Applications</b> | ELISA, IF; Recommended dilution: IF:1:50-1:200                                          |
| <b>Form</b>                | Liquid                                                                                  |
| <b>Conjugate</b>           | Non-conjugated                                                                          |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4        |
| <b>Purification Method</b> | >95%, Protein G purified                                                                |
| <b>Isotype</b>             | IgG                                                                                     |
| <b>Clonality</b>           | Polyclonal                                                                              |
| <b>Alias</b>               | T-box transcription factor TBX1 (T-box protein 1) (Testis-specific T-box protein), TBX1 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                    |
| <b>Research Area</b>       | Neuroscience                                                                            |
| <b>Target Names</b>        | TBX1                                                                                    |

## Image



Immunofluorescent analysis of PC-3 cells using CSB-PA023245YA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)