



CD36 Antibody, HRP conjugated

Product Code	CSB-PA09479B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16671
Immunogen	Recombinant Human Platelet glycoprotein 4 protein (191-376AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Platelet glycoprotein 4 (Fatty acid translocase) (FAT) (Glycoprotein IIIb) (GPIIIB) (Leukocyte differentiation antigen CD36) (PAS IV) (PAS-4) (Platelet collagen receptor) (Platelet glycoprotein IV) (GPIV) (Thrombospondin receptor) (CD antigen CD36), CD36, GP3B GP4
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CD36
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.