



TGFB111 Antibody, FITC conjugated

Product Code	CSB-PA13049C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43294
Immunogen	Recombinant Human Transforming growth factor beta-1-induced transcript 1 protein (138-200AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming growth factor beta-1-induced transcript 1 protein (Androgen receptor coactivator 55 kDa protein) (Androgen receptor-associated protein of 55 kDa) (Hydrogen peroxide-inducible clone 5 protein) (Hic-5), TGFB111, ARA55
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TGFB111
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.