



AKT3 Antibody, FITC conjugated

Product Code	CSB-PA15929C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y243
Immunogen	Recombinant Human RAC-gamma serine/threonine-protein kinase protein (1-164AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RAC-gamma serine/threonine-protein kinase (EC 2.7.11.1) (Protein kinase Akt-3) (Protein kinase B gamma) (PKB gamma) (RAC-PK-gamma) (STK-2), AKT3, PKBG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AKT3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.