



PAX5 Antibody, Biotin conjugated

Product Code	CSB-PA20659D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02548
Immunogen	Recombinant Human Paired box protein Pax-5 protein (362-391AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Paired box protein Pax-5 (B-cell-specific transcription factor) (BSAP), PAX5
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PAX5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.