



HLA-C Antibody, HRP conjugated

Product Code	CSB-PA320739LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10321
Immunogen	Recombinant Human HLA class I histocompatibility antigen, Cw-7 alpha chain protein (62-306AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class I histocompatibility antigen, Cw-7 alpha chain (MHC class I antigen Cw*7), HLA-C, HLAC
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HLA-C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.