



WDPCP Antibody, HRP conjugated

Product Code	CSB-PA528516LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95876
Immunogen	Recombinant Human WD repeat-containing and planar cell polarity effector protein fritz homolog protein (621-737AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat-containing and planar cell polarity effector protein fritz homolog (hFRTZ) (Bardet-Biedl syndrome 15 protein) (WD repeat-containing and planar cell polarity effector protein), WDPCP, BBS15 C2orf86 FRITZ
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WDPCP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.