



WDPCP Antibody, Biotin conjugated

Product Code	CSB-PA528516LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95876
Immunogen	Recombinant Human WD repeat-containing and planar cell polarity effector protein fritz homolog protein (621-737AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat-containing and planar cell polarity effector protein fritz homolog (hFRTZ) (Bardet-Biedl syndrome 15 protein) (WD repeat-containing and planar cell polarity effector protein), WDPCP, BBS15 C2orf86 FRITZ
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WDPCP