




HLA-B Antibody

Product Code	CSB-PA642773LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q29718
Immunogen	Recombinant Human HLA class I histocompatibility antigen, B-82 alpha chain protein (197-306AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class I histocompatibility antigen, B-82 alpha chain (MHC class I antigen B*82), HLA-B, HLAB
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HLA-B

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate, MCF-7 whole cell lysate


All lanes: HLA-B antibody at 3.2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 41 kDa

Observed band size: 41 kDa



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA642773LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.