



HLA-B Antibody, Biotin conjugated

Product Code	CSB-PA642773LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q29718
Immunogen	Recombinant Human HLA class I histocompatibility antigen, B-82 alpha chain protein (197-306AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class I histocompatibility antigen, B-82 alpha chain (MHC class I antigen B*82), HLA-B, HLAB
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HLA-B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.