



WWC2 Antibody, HRP conjugated

Product Code	CSB-PA734879LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6AWC2
Immunogen	Recombinant Human Protein WWC2 protein (828-1127AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein WWC2 (BH-3-only member B) (WW domain-containing protein 2), WWC2, BOMB
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WWC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.