



ASPM Antibody, Biotin conjugated

Product Code	CSB-PA814245LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZT6
Immunogen	Recombinant Human Abnormal spindle-like microcephaly-associated protein (101-234AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Abnormal spindle-like microcephaly-associated protein (Abnormal spindle protein homolog) (Asp homolog), ASPM, MCPH5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ASPM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.