



MCAT Antibody, FITC conjugated

Product Code	CSB-PA815557LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IVS2
Immunogen	Recombinant Human Malonyl-CoA-acyl carrier protein transacylase, mitochondrial protein (22-117AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Malonyl-CoA-acyl carrier protein transacylase, mitochondrial (MCT) (EC 2.3.1.39) (Mitochondrial malonyl CoA:ACP acyltransferase) (Mitochondrial malonyltransferase) ([Acyl-carrier-protein] malonyltransferase), MCAT, MT
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MCAT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.