



HSCB Antibody, HRP conjugated

Product Code	CSB-PA816893LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWL3
Immunogen	Recombinant Human Iron-sulfur cluster co-chaperone protein HscB, mitochondrial protein (30-235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Iron-sulfur cluster co-chaperone protein HscB, mitochondrial (DnaJ homolog subfamily C member 20) (Hsc20), HSCB, DNAJC20 HSC20
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSCB