



HSCB Antibody, FITC conjugated

Product Code	CSB-PA816893LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWL3
Immunogen	Recombinant Human Iron-sulfur cluster co-chaperone protein HscB, mitochondrial protein (30-235AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Iron-sulfur cluster co-chaperone protein HscB, mitochondrial (DnaJ homolog subfamily C member 20) (Hsc20), HSCB, DNAJC20 HSC20
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSCB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.