



NIN Antibody, HRP conjugated

Product Code	CSB-PA843155LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N4C6
Immunogen	Recombinant Human Ninein protein (80-185AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ninein (hNinein) (Glycogen synthase kinase 3 beta-interacting protein) (GSK3B-interacting protein), NIN, KIAA1565
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NIN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.