



ORAI1 Antibody, Biotin conjugated

Product Code	CSB-PA846601LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96D31
Immunogen	Recombinant Human Calcium release-activated calcium channel protein 1 protein (1-87AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium release-activated calcium channel protein 1 (Protein orai-1) (Transmembrane protein 142A), ORAI1, CRACM1 TMEM142A
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ORAI1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.