



MCU Antibody, FITC conjugated

Product Code	CSB-PA847694LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NE86
Immunogen	Recombinant Human Calcium uniporter protein, mitochondrial protein (155-233AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium uniporter protein, mitochondrial (HsMCU) (Coiled-coil domain-containing protein 109A), MCU, C10orf42 CCDC109A
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MCU