



SRMS Antibody, Biotin conjugated

Product Code	CSB-PA862011LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H3Y6
Immunogen	Recombinant Human Tyrosine-protein kinase Srms protein (19-155AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein kinase Srms (EC 2.7.10.2), SRMS, C20orf148
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SRMS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.