




DEFB126 Antibody

Product Code	CSB-PA866323LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYW3
Immunogen	Recombinant Human Beta-defensin 126 protein (36-111AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-defensin 126 (Beta-defensin 26) (DEFB-26) (Defensin, beta 126) (Epididymal secretory protein 13.2) (ESP13.2) (HBD26), DEFB126, C20orf8 DEFB26
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DEFB126

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA866323LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.