



DEFB126 Antibody, Biotin conjugated

Product Code	CSB-PA866323LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYW3
Immunogen	Recombinant Human Beta-defensin 126 protein (36-111AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-defensin 126 (Beta-defensin 26) (DEFB-26) (Defensin, beta 126) (Epididymal secretory protein 13.2) (ESP13.2) (HBD26), DEFB126, C20orf8 DEFB26
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DEFB126
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.