




MENT Antibody

Product Code	CSB-PA871579LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BUN1
Immunogen	Recombinant Human Protein MENT protein (24-197AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein MENT (Methylated in normal thymocytes protein), MENT, C1orf56
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MENT

Image



Western Blot

Positive WB detected in: MCF-7 whole cell lysate, Rat brain tissue


All lanes: MENT antibody at 3.2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 37, 19 kDa

Observed band size: 37 kDa



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA871579LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.