



MENT Antibody, HRP conjugated

Product Code	CSB-PA871579LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BUN1
Immunogen	Recombinant Human Protein MENT protein (24-197AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein MENT (Methylated in normal thymocytes protein), MENT, C1orf56
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MENT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.