



PDF Antibody, Biotin conjugated

Product Code	CSB-PA872538LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBH1
Immunogen	Recombinant Human Peptide deformylase, mitochondrial protein (16-143AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptide deformylase, mitochondrial (EC 3.5.1.88) (Polypeptide deformylase), PDF, PDF1A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.