



PSORS1C1 Antibody, HRP conjugated

Product Code	CSB-PA883431LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIG5
Immunogen	Recombinant Human Psoriasis susceptibility 1 candidate gene 1 protein (1-152AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Psoriasis susceptibility 1 candidate gene 1 protein (Protein SEEK1), PSORS1C1, C6orf16 SEEK1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSORS1C1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.