



PSORS1C1 Antibody, Biotin conjugated

Product Code	CSB-PA883431LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIG5
Immunogen	Recombinant Human Psoriasis susceptibility 1 candidate gene 1 protein (1-152AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Psoriasis susceptibility 1 candidate gene 1 protein (Protein SEEK1), PSORS1C1, C6orf16 SEEK1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSORS1C1