



POLE Antibody, Biotin conjugated

Product Code	CSB-PA894382LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07864
Immunogen	Recombinant Human DNA polymerase epsilon catalytic subunit A protein (1175-1311AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA polymerase epsilon catalytic subunit A (EC 2.7.7.7) (DNA polymerase II subunit A), POLE, POLE1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POLE