



POLE Antibody, Biotin conjugated

Product Code	CSB-PA894382LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07864
Immunogen	Recombinant Human DNA polymerase epsilon catalytic subunit A protein (1175-1311AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA polymerase epsilon catalytic subunit A (EC 2.7.7.7) (DNA polymerase II subunit A), POLE, POLE1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POLE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.