



WNK2 Antibody, HRP conjugated

Product Code	CSB-PA896509LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3S1
Immunogen	Recombinant Human Serine/threonine-protein kinase WNK2 protein (1177-1464AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase WNK2 (EC 2.7.11.1) (Antigen NY-CO-43) (Protein kinase lysine-deficient 2) (Protein kinase with no lysine 2) (Serologically defined colon cancer antigen 43), WNK2, KIAA1760 PRKWINK2 SDCCAG43
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WNK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.