



MYRF Antibody, Biotin conjugated

Product Code	CSB-PA897527LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2G1
Immunogen	Recombinant Human Myelin regulatory factor protein (843-1004AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myelin regulatory factor (EC 3.4.-.-) (Myelin gene regulatory factor) [Cleaved into: Myelin regulatory factor, N-terminal; Myelin regulatory factor, C-terminal], MYRF, C11orf9 KIAA0954 MRF
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MYRF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.