




AURKC Antibody


Product Code	CSB-PA002460LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQB9
Immunogen	Recombinant Human Aurora kinase C protein (12-143AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aurora kinase C (EC 2.7.11.1) (Aurora 3) (Aurora/IPL1-related kinase 3) (ARK-3) (Aurora-related kinase 3) (Aurora/IPL1/Eg2 protein 2) (Serine/threonine-protein kinase 13) (Serine/threonine-protein kinase aurora-C), AURKC, AIE2 AIK3 AIRK3 ARK3 STK13
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	AURKC

Image




Western Blot


Positive WB detected in: A375 whole cell lysate
 All lanes: AURKC antibody at 3.2µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 36, 33, 34 kDa
 Observed band size: 45 kDa



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA002460LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human melanoma using CSB-PA002460LA01HU at dilution of 1:100



Immunofluorescent analysis of Hela cells using CSB-PA002460LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.