



notch1a Antibody, Biotin conjugated

Product Code	CSB-PA015949LD01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46530
Immunogen	Recombinant Zebrafish Neurogenic locus notch homolog protein 1 protein (21-338AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurogenic locus notch homolog protein 1 (Notch 1) [Cleaved into: Notch 1 extracellular truncation (NEXT); Notch 1 intracellular domain (NICD)], notch1a, notch
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Neuroscience
Target Names	notch1a
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.