



tbx5a Antibody, FITC conjugated

Product Code	CSB-PA881254LC01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9IAK8
Immunogen	Recombinant Zebrafish T-box transcription factor TBX5-A protein (1-492AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-box transcription factor TBX5-A (zftbx5a) (T-box protein 5) (zTbx5), tbx5a, heartstrings hst tbx5 tbx5.1
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Epigenetics and Nuclear Signaling
Target Names	tbx5a
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.