



Recombinant Yellow fever virus Genome polyprotein, partial

Product Code	CSB-MP365905YAC1
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P03314
Product Type	Recombinant Protein
Immunogen Species	Yellow fever virus (strain 17D vaccine) (YFV)
Purity	>85% (SDS-PAGE)
Sequence	AHCIGITDRDFIEGVHGGTWVSATLEQDKCVTVMAPDKPSLDISLETVAIDRPA EVRKVCYNAVLTHVKINDKCPSTGEAHLAEENEGDNACKRTYSDRGWGNGC GLFGKGSIVACAKFTCAKSMSLFEVDQTKIQYVIRAQLHVGAKQENWNTDIKTL KFDALSGSQEVEFIGYGKATLECQVQTAVDGNSYIAEMETESWIVDRQWAQ DLTLPWQSGSGGVWREMHHLVEFEPHAATIRVLALGNQEGSLKTALTGAMR VTKDTNDNLYKLHGGHVSCRVKLSALTLKGTSYKICTDKMFFVKNPTDTGHG TVVMQVKVSKGAPCRIPVIVADDLTAANKGILVTVNPIASTNDDEVLIENVPPFG DSYIIVGRGDSRLTYQWHKEGSSIGKLFTQTMKGVERLAVMGDTAWDFSSAG GFFTSVGKGIHTVFGSAFQGL
Source	Mammalian cell
Protein Names	Recommended name: Genome polyprotein Cleaved into the following 14 chains: 1. Capsid protein C Alternative name(s): Core protein prM Peptide pr Small envelope protein M Alternative name(s): Matrix protein Envelope protein E Non-str
Expression Region	286-730aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.