




Recombinant Plasmodium berghei Major facilitator superfamily domain-containing protein,partial

Product Code	CSB-CF2221PLHa2
Abbreviation	Recombinant Plasmodium berghei PBK173_000211300 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A0A0Y9WQZ1
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose
Product Type	Transmembrane Proteins
Immunogen Species	Plasmodium berghei
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	QNIKTIFSSFTFYFFFFFIFIYAFLLSIMHIFINYFFYIYLFVFNINIIYISNIYTIIMTFASLI AIPFSGYIIDNIGSFLFLLLCSSFFILIAISGTIYSCVFNLRSEVIAFISFNLIGISESIIP TVIISQIPTHLCVKKNEDITAAFAIFELVSMLIVSVNNYIFGYFLINKEYLNGLYILFV FVILVISLIFLLIFTIYWKAR
Research Area	others
Source	in vitro E.coli expression system
Target Names	N/A
Expression Region	756-960aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	39.9kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.