




Recombinant Human Cytochrome b-245 heavy chain (CYBB), partial

Product Code	CSB-EP006325HU1a2
Relevance	Critical component of the mbrane-bound oxidase of phagocytes that generates superoxide. It is the terminal component of a respiratory chain that transfers single electrons from Cytoplasmic domain NADPH across the plasma mbrane to molecular oxygen on the exterior. Also functions as a voltage-gated proton channel that mediates the H ⁺ currents of resting phagocytes. It participates in the regulation of cellular pH and is blocked by zinc.
Abbreviation	Recombinant Human CYBB protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04839
Alias	CGD91-phoxCytochrome b(558) subunit beta ;Cytochrome b558 subunit betaHeme-binding membrane glycoprotein gp91phoxNADPH oxidase 2Neutrophil cytochrome b 91 kDa polypeptide;Superoxide-generating NADPH oxidase heavy chain subunitgp91-1gp91-phoxp22 phagocyte B-cytochrome
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	ERLVRFWRSQQKVVITKVVTHPFTKIELQMKKKGFKMEVGQYIFVKCPKVS KL EWHPFTLTSAPEEDFFSIHIRIVGDWTEGLFNACGCDKQEFQDAWKLPKIAVD GPFGTASEDVFSYEVVMLVGAGIGVTPFASILKSVWYKYCNNATNLKLKKIYFY WLCRDTHAFEWFADLLQLLESQMQRNNAGFLSYNIYLTGWDESQANHFAVH HDEEKDVITGLKQKTLYGRPNWDNEFKTIASQHPNTRIGVFLCGPEALAE TL SK QSI SNSESGPRGVHFIFNKENF
Research Area	Cancer
Source	E.coli
Target Names	CYBB
Expression Region	283-570aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	49.2kDa
Protein Length	Cytoplasmic Domain



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.