




Recombinant Ceratopteris richardii Cyanovirin-N homolog

Product Code	CSB-YP309974CGPe1
Relevance	Mannose-binding lectin
Abbreviation	Recombinant Ceratopteris richardii Cyanovirin-N homolog protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P86326
Product Type	Recombinant Proteins
Immunogen Species	Ceratopteris richardii (Triangle waterfern)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	QCNFANSCTGVELYGYILRGDCINEDGHPHATSINLNYYIGNDNGRLEYPGES FGSSCVKTALNDGHTLTASCKGADGQYHDSSMDLNYVVGNSYGYMEPCRAS NADHVLKSSSE
Research Area	Others
Source	Yeast
Target Names	N/A
Expression Region	28-142aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	12.4kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.