




# Recombinant Human Sulfotransferase 1A3 (SULT1A3)

<b>Product Code</b>	CSB-EP022935HUe0
<b>Relevance</b>	Sulfotransferase that utilizes 3'-phospho-5'-adenylyl sulfate (PAPS) as sulfonate donor to catalyze the sulfate conjugation of phenolic monoamines (neurotransmitters such as dopamine, norepinephrine and serotonin) and phenolic and catechol drugs.
<b>Abbreviation</b>	Recombinant Human SULT1A3 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0DMM9
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	MELIQDTSRPPLEYVKGVPLIKYFAEALGPLQSFQARPDDLINTYPKSGTTWV SQILDMIYQGGDLEKCNRAPIYVRVPFLEVNDPGEPSGLETCLKDTPPPRLIKSH LPLALLPQTLDDQKVKVYVARNPKDVAVSYYHFHRMEKAHPEPGTWDSFLE KFMAGEVSYGSWYQHVVQEWELSRTHPVLVLYFYEDMKENPKREIQKILEFVG RSLPEETMDFMVQHTSFKEMKKNPMTNYTTVPQELMDHSISPFRKGMAGD WKTFTTVAQNERFDADYAEKMAGCSLSFRSEL
<b>Research Area</b>	Neuroscience
<b>Source</b>	E.coli
<b>Target Names</b>	SULT1A3
<b>Expression Region</b>	1-295aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal GST-tagged
<b>Mol. Weight</b>	61.2kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.