




Recombinant Juniperus ashei Pectate lyase 1

Product Code	CSB-YP305864JANe1
Relevance	Has low pectate lyase activity.
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P81294
Alias	Major pollen allergen Jun a 1 Allergen: Jun a 1
Product Type	Recombinant Protein
Immunogen Species	Juniperus ashei (Ozark white cedar)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	DNPIDSCWRGDSNWDQNRMKLADCAVGFSGSSTMGGKGGDFYTVTSTDDNP VNPTPGTLRYGATREKALWIIFSQNMNIKLMPLYVAGHKTIDGRGADVHLGN GGPCLFMRKVSHVILHSLHIHGCNTSVLGDVLVSESIGVEPVHAQDGDAITMR NVTNAWIDHNSLSDCSDGLIDVTLGSTGITISNNHFFNHHKVMMLLGHDDTYDDD KSMKVTVAFNQFGPNAGQRMPRARYGLVHVANNNYDPWNIYAIGSSNPTIL SEGNSFTAPSESYKKEVTKRIGCESPSACANWVWRSTRDAFINGAYFVSSGK TEETNIYNSNEAFKVENGNAAPQLTKNAGVVT
Research Area	Others
Source	Yeast
Expression Region	22-367aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	37.6kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.